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
Monaco 5.51 Online Help

Contents Index Search

- What are the Monaco Activities?
- What is New in this Monaco Release?
- Copyright and Trademarks
- Hardware and Software Limitations
- Glossary of Terms
- Getting Started
- Ribbons
- Layouts
- Opening a Patient
- Workspace Control
- Remove Patient Dialog Box
- Keyboard Shortcuts
- About the Tool Controls
- Why is a menu command or button visible?
- Positioning Screen Elements
- Using Elekta Resources
- Downloading Software Patches
- About Your Customized Settings
- Reset Control Positions
- Log
- Patient Access Log
- References
- Playing Flash Movies
- Exit from Monaco
- Image Manipulation
- Tools
- Planning Activity
- Monaco Extension Architecture

References

Monaco General References

Alber, M PhD-thesis, University of Tübingen, <http://w210.ub.uni-tuebingen.de/dbt/volltexte/2001/221>.

Alber M, Paulsen F, Eschmann SM, Machulla HJ. On biologically conformal boost dose optimization. Phys Med Biol. 2005; 50(12): 3229-39. PMID: 12587912 [PubMed - indexed for MEDLINE]

Alber M, Meedt G, Nusslin F, Reemtsen R. On the degeneracy of the IMRT optimization problem. Med Phys. 2005; 32(12): 3229-39. PMID: 12462725 [PubMed - indexed for MEDLINE]

Alber M, Birkner M, Nusslin F. Tools for the analysis of dose optimization: II. Sensitivity analysis. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 12408484 [PubMed - indexed for MEDLINE]

Alber M, Nusslin F. Tools for the analysis of dose optimization: I. Effect-volume histogram. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 12171333 [PubMed - indexed for MEDLINE]

Alber M, Nusslin F. [A concept for the optimization of clinical IMRT] Z Med Phys. 2002; 12(2): 109-13. G-indexed for MEDLINE]

Alber M, Nusslin F. Optimization of intensity modulated radiotherapy under constraints for static and dynamic IMRT. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 11768502 [PubMed - indexed for MEDLINE]


Alber M, Nusslin F. A representation of an NTCP function for local complication mechanisms. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 11229724 [PubMed - indexed for MEDLINE]

Alber M, Nusslin F. Intensity modulated photon beams subject to a minimal surface smoothing constraint. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 10843115 [PubMed - indexed for MEDLINE]

Alber M, Nusslin F. On the appropriateness of the probability model for series type complications. Radiother Oncol. 2006; 80(2): 109-13. PMID: 10577691 [PubMed - indexed for MEDLINE]

Alber M, Nusslin F. An objective function for radiation treatment optimization based on local biological models. Phys Med Biol. 2006; 51(12): 3229-39. PMID: 10070796 [PubMed - indexed for MEDLINE]

Patient Access Log

Open Patient

Close Patient

Delete Patients

Save

Save Plan As

Save Template

Save Template As

Save Summation/Difference Plan

Manage Templates

Settings

Setup

Treatment Couch

Log Files

Monaco application and patient transfer Logfiles

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Access Log Viewer [PROSTATE Fusion Prostate]


File Edit View

☒ Auto Update ☒ Auto Scroll


Time	(All)	User Name	(All)	Event
4	12/2/2019 10:37:59 AM	FOCUS		Loaded:: Patient.
5	12/2/2019 10:38:08 AM	FOCUS		Loaded:: Studysset: CT1; StructureSet: SS_CT1.
6	12/6/2019 10:22:15 AM	FOCUS		Loaded:: Patient.
7	12/6/2019 10:22:25 AM	FOCUS		Loaded:: Studysset: CTClean; StructureSet: SS_CTClean.
8	12/6/2019 10:23:51 AM	FOCUS		Saved:: Studysset: CTClean; StructureSet: SS_CTClean.
9	12/6/2019 10:25:19 AM	FOCUS		Saved plan as:: Studysset: CTClean; StructureSet: SS_CTClean; Plan: 1.
10	12/6/2019 10:25:25 AM	FOCUS		Deleted:: Studysset: CTClean; StructureSet: SS_CTClean; Plan: 1.
11	12/6/2019 10:25:28 AM	FOCUS		Loaded:: Studysset: CTClean; StructureSet: SS_CTClean.
12	12/6/2019 12:01:45 PM	FOCUS		Loaded:: Patient.

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
All Activities

Save the patient data	Ctrl + S
Undo the last contour or block editing action	Ctrl + Z
Cut the selected contour	Ctrl + X
Copy the selected contour	Ctrl + C
Paste the cut or copied contour	Ctrl + V
Exit Monaco	Ctrl + Q or Alt + F4
Print BEV DRR	Ctrl + P
Open online help	F1
View next superior/inferior image	Page Up/Page Down
Zoom in	Numeric keypad +
Zoom out	Numeric keypad -
Print entire screen	Ctrl + Prnt Scrn
Print active window only into the clipboard	Ctrl + Alt + PrntScrn

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Planning Activity	
DICOM Export	Ctrl + E
Start a new plan	Ctrl + N
Close a plan	Ctrl + W
Delete a plan	Ctrl + D
Open a patient	Ctrl + O
Copy selected contour to the superior/inferior image	Ctrl + Page Up (copy superior) Ctrl + Page Down (copy inferior)
Delete the last drawn contour point or section	Backspace
Delete a selected contour or port object	Delete
Remove the current selected structure from the studyset	Shift + Delete
Quit an incomplete drawing or editing operation	Esc
Increase guide radius value	→ or ↑ or >
Decrease guide radius value	← or ↓ or <

Room's Eye View	
Toggle full screen mode	F
Show couchtop	U
Hide couchtop	I
Toggle drawing couchtop as lines	O
Toggle drawing couchtop as points	P
Toggle gantry visibility	G
Toggle décor visibility	D
Toggle room lasers on/off	L
Toggle axes on/off	A
Command	Shortcut
Scroll segments in Fluence View	Ctrl + left/right arrow

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
User Name : **elekta**
Password : **focus1.1**

Monaco

A screenshot of a "User Validation" dialog box. It has a dark blue header. Inside, there are two input fields: "User Name:" with the value "focus" and "Password:" with masked characters "*****". Below the fields are two buttons: "OK" and "Cancel".

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DCM	Patient Name	Patient ID	Creation Date
4DLung, TRN	4DLung		Monday, April 29, 2013
FeetFirst, Supine	FeetFirstLeg		Thursday, December 27, 2012
SUVCalculation, TRN	CTPETLung		Monday, October 15, 2012
Demo, Patient	DynConform		Saturday, September 11, 2010
Prostate, Nodes	ProstateNodes		Thursday, August 26, 2010
Bias, Dose	WeightLossTRN		Thursday, August 5, 2010
Brain, STargets	BrainSTrgtsTRN		Sunday, June 20, 2010
LeftLung, SBRT	LLungSBRTTRN		Friday, June 4, 2010
Fusion, Prostate	PROSTATE		Wednesday, August 19, 2009
Patient, HIN	HeadandNeck		Tuesday, May 26, 2009
ESOPHAGUS, TRN	Esophagus		Tuesday, January 2, 2007
HINPlan, TRN	HINPlan		Tuesday, November 8, 2005
AdvHINTRN, Head	AdvancedHIN		Friday, June 24, 2005
Monaco, Training	TestPatient		

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
Structures												
View: Contoured All Layers Adapt S												
Name	Color	Visible	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlin...	2D Transparency	3D/BEV Transparency		
BLADDER		<input checked="" type="checkbox"/>	352.494	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
CTV		<input checked="" type="checkbox"/>	57.345	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
GTV		<input checked="" type="checkbox"/>	43.011	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
patient		<input checked="" type="checkbox"/>	13757.324	External	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
PTV		<input checked="" type="checkbox"/>	122.356	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
RECTUM		<input checked="" type="checkbox"/>	63.420	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
SV		<input checked="" type="checkbox"/>	11.281	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>				
<click to add a new row>												

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
Structures

View: Contours All Layers Adapt Setup

Name	Color	Visible	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlin...	2D Transparency	3D/BEV Transparency
BLADDER	Yellow	<input checked="" type="checkbox"/>	352.494	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CTV	Blue	<input checked="" type="checkbox"/>	57.345	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Range: 0.010 to 15.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
GTV	Red	<input checked="" type="checkbox"/>	43.011	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
patient	Brown	<input checked="" type="checkbox"/>	13757.324	External	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
PTV	Purple	<input checked="" type="checkbox"/>	122.356	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
RECTUM	Green	<input checked="" type="checkbox"/>	63.420	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SV	Pink	<input checked="" type="checkbox"/>	11.281	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Range: 0.010 to 15.000

Monte Carlo

Collapsed Cone and eMC

Structures


View: Contours All Layers Adapt Setup

Name	Color	Visible	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlin...	2D Transparency	3D/BEV Transparency
BLADDER	Yellow	<input checked="" type="checkbox"/>	352.494	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CTV	Blue	<input checked="" type="checkbox"/>	57.345	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Range: 0.010 to 2.456	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
GTV	Red	<input checked="" type="checkbox"/>	43.011	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
patient	Brown	<input checked="" type="checkbox"/>	13757.324	External	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
PTV	Purple	<input checked="" type="checkbox"/>	122.356	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
RECTUM	Green	<input checked="" type="checkbox"/>	63.420	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SV	Pink	<input checked="" type="checkbox"/>	11.281	Internal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

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


Anatomical Groupの作成

ToolsタブのAnatomical Groupsをクリック

Anatomical Groups													
Anatomical Group Name: Chest													
Name	Color	Type	Force ED on CT	Force ED on MR	Force Relative ED	Fill ED on CT	Fill Relative ED	Layer	Show 2D Outlines	2D Transparency	3D/BEV Transparency	Generation	
LEFT LUNG		Internal	<input type="checkbox"/>	<input checked="" type="checkbox"/>	1.000	<input type="checkbox"/>		3	<input checked="" type="checkbox"/>			Deform	
RIGHT LUNG		Internal	<input type="checkbox"/>	<input checked="" type="checkbox"/>	1.000	<input type="checkbox"/>		1	<input checked="" type="checkbox"/>			Deform	
SPINAL CORD		Internal	<input type="checkbox"/>	<input checked="" type="checkbox"/>	1.000	<input type="checkbox"/>		2	<input checked="" type="checkbox"/>			Rigid	
<click to add a new row>													
<div> <div>Create New Anatomical Group</div> <div>Copy Anatomical Group</div> <div>Delete Anatomical Group</div> <div>Apply</div> <div>OK</div> <div>Cancel</div> </div>													

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Structures		
View:	Contoured	All Layers Adapt Setup
Name	Color	Generation
LT LUNG	Blue	Deform
CORD	Blue	Deform
ESOPHAGUS	Green	Deform
HEART	Red	Rigid
LIVER	Orange	User
LT KIDNEY	Dark Blue	Deform
RT KIDNEY	Light Green	Deform
RT LUNG	Yellow	Deform
BOWEL	Dark Green	Deform

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Structures						
View:	Contoured	All	Layers	Adapt Setup		
Name	Color	Generation	Force ED on MR	Relative ED	Mean ED	
PTV	Blue	Margin	<input checked="" type="checkbox"/>	1.000	1.000	
CTV	Red	Margin	<input checked="" type="checkbox"/>	1.000	0.981	

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Structures								
View:	Contoured	All	Layers	Adapt Setup				
Name	Color	Force ED	Fill ED	Relative ED				
GTV	Red	<input type="checkbox"/>	<input type="checkbox"/>					
SV	Purple	<input type="checkbox"/>	<input type="checkbox"/>					
BLADDER	Yellow	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.000				
RECTUM	Green	<input type="checkbox"/>	<input type="checkbox"/>					
CTV	Cyan	<input type="checkbox"/>	<input type="checkbox"/>					
PTV	Blue	<input type="checkbox"/>	<input type="checkbox"/>					
patient	Brown	<input type="checkbox"/>	<input type="checkbox"/>					

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Structures								
View:	Contoured	All	Layers	Adapt Setup				
Name	Color	Visible	Volume (cm ³)	Type	Generation	Force ED	Fill ED	Relative ED
BLADDER	Yellow	<input checked="" type="checkbox"/>	352.232	Internal	Deform	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.000
CTV	Cyan	<input checked="" type="checkbox"/>	54.145	Internal	Margin	<input type="checkbox"/>	<input type="checkbox"/>	
GTV	Red	<input checked="" type="checkbox"/>	42.896	Internal	Deform	<input type="checkbox"/>	<input type="checkbox"/>	
patient	Brown	<input checked="" type="checkbox"/>	13640.678	External	Deform	<input type="checkbox"/>	<input type="checkbox"/>	
PTV	Blue	<input checked="" type="checkbox"/>	121.647	Internal	Margin	<input type="checkbox"/>	<input type="checkbox"/>	
RECTUM	Green	<input checked="" type="checkbox"/>	63.149	Internal	Deform	<input type="checkbox"/>	<input type="checkbox"/>	
SV	Purple	<input checked="" type="checkbox"/>	11.250	Internal	Deform	<input type="checkbox"/>	<input type="checkbox"/>	

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Structures								
View:	Contoured	All	Layers	Adapt Setup				
Name	Color	Generation	Force ED on MR	Relative ED	←	Mean ED		
GTV	Red	Deform	<input checked="" type="checkbox"/>	1.072		1.072		
SV	Purple	Deform	<input checked="" type="checkbox"/>	1.062		1.062		
BLADDER	Yellow	Deform	<input checked="" type="checkbox"/>	1.000	<input type="text" value="1.000"/>	1.202		
RECTUM	Green	Deform	<input checked="" type="checkbox"/>	0.960		0.960		
CTV	Cyan	Margin	<input checked="" type="checkbox"/>	1.070		1.070		
PTV	Blue	Margin	<input checked="" type="checkbox"/>	1.062		1.062		
patient	Brown	Deform	<input checked="" type="checkbox"/>	1.047		1.047		

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
Treatment Unit not approved for Clinical use: Testing
Electron densities are overridden on structures that may be overlapped
The displayed electron density values are for adapt setup visualization only

Name	Color	Force ED	Fill ED	Relative ED
GTV		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.072
SV		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.062
CTV		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.070
RECTUM		<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.980
PTV		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.062
BLADDER		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.000
patient		<input checked="" type="checkbox"/>	<input type="checkbox"/>	1.047

Grid Type: Adapt Setup Electr Structure:
Volume Cursor: SS Units: 0.900 1.100
Grid Volume

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<click to add a new row>

Create New Anatomical Group Copy Anatomical Group Delete Anatomical Group Apply OK Cancel

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Name: NewTempPlan Description:

Delivery: VMAT Anatomical Site: All

Scan Orientation (ANUS): Head First Supine

Treatment Orientation: ☒ Head First ☐ Feet First

Select template to import

- ☐ Template: DEFAULTVMAT (Rx Site: , Rx Dose: cGy, Total Beams: 1)
- ☐ Template: RTOG0529AnalCanal (Rx Site: 1, Rx Dose: 5040.0 cGy, Total Beams: 1)
- ☐ Template: RTOG0615NPC (Rx Site: , Rx Dose: 7000.0 cGy, Total Beams: 1)
- ☒ Template: RTOG0631SpineVMAT (Rx Site: Spine, Rx Dose: 1600.0 cGy, Total Beams: 1)
- ☒ Vmat (Number of Arcs: 1)
- ☐ Template: RTOG0813LUNG (Rx Site: Lung, Rx Dose: 5000.0 cGy, Total Beams: 1)
- ☐ Template: RTOG0815Prostate (Rx Site: , Rx Dose: 7920.0 cGy, Total Beams: 1)

MOSAIQ Options

Course ID: 1 Plan Intent: Curative Tolerance Table: 1 Clinical

Beam	Treatment Unit	Map Machine	Modality	Algorithm	Energy	Isocenter Location	X(cm)	Y(cm)	Z(cm)
1	agility		Photon	Monte Carlo	6.0 MV	Center of PTV_4500	-0.73	-1.55	0.92

Port Options

☒ Import Beams Only ☐ Retain Template Beam Shapes ☐ Auto-conform Ports

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Name: Plan1 Description:

Delivery: 3D Anatomical Site: All

Scan Orientation (CT1): Head First Supine

Treatment Orientation

Head First Feet First

MONACO Options

Monaco

The selected treatment orientation (Feet First) does not match the orientation in the images (Head First). Do you want to proceed?

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
※Dosimetric Criteria, Reference Dose, IMRT Constraintなどに適用

Map Structures

Some structures that are referenced in the template are not available in the selected structure set. Dosimetric criteria and reference doses relating to structures that are not mapped will be discarded.

Template structure	Available structure
PTV	PTV66

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Beams

Delete Parent Beams


Beam	Description	Field ID	Visible	Delivery	Treatment Unit	Modality	Algorithm	Energy
1	AP	1	✓	3D	Elekta	Photon	Collapsed Cone	6.0 MV
2	New beam	2	✓	3D	Elekta	Photon	Collapsed Cone	6.0 MV
3	New beam	3	✓	3D	Elekta	Photon	Collapsed Cone	6.0 MV

<click to add a new beam>

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Structures Prescription **Beams** Dose Reference Points DVH Statistics

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Current Port: 0

Port: [Block] [MLC]

Auto Conform
Structure: PTV
Margin (cm): 1.00

Apply to Visible Beams

☐ Snap Jaws to Port

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Auto Conform
Structure: PTV
Margin (cm): 2.00

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
☒ Snap Jaws to Port
☐ Rotate Port with Collimator
 Port by Shapes: NONE
 Editing Tool: ☒ Replace ☐ Reshape

Editing Tool: ☒ Replace ☐ Reshape

Pair #	Left	Right	Parked
38	-5.00	5.00	<input type="checkbox"/>
39	-5.00	2.00	<input type="checkbox"/>
40	-5.00	2.00	<input type="checkbox"/>
41	-5.00	2.00	<input type="checkbox"/>
42	-5.00	2.00	<input type="checkbox"/>
43	-5.00	5.00	<input type="checkbox"/>

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Prescription

Prescription Segments

Add Rx Delete Rx

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Rx ID Rx Site Prescribe To Rx Dose (cGy) Number of Fractions Fractional Dose (cGy)

▼ Physician's Intent A ▼ Plan Isocenter X 0.50 Y -121.70 Z -0.33 7000.0 35 200.0

Actual Dose = 0.0 cGy

Rescale 7000.0 cGy to...

Weight beams by: ☒ Dose ☐ MU

均等配分 Equal Weights

Beam	Description	Field ID	%	Lock	MU / Fx
1 g0		1	25.00	<input type="checkbox"/>	0.00
2 g90		2	25.00	<input type="checkbox"/>	0.00
3 g180		3	25.00	<input type="checkbox"/>	0.00
4 g270		4	25.00	<input type="checkbox"/>	0.00

Total MU / Fx 0.00


Structure: Prescription Beams Dose Reference Points DVH Statistics

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The screenshot shows the 'Prescription' window in the Elekta software. A dropdown menu is open, listing various anatomical points and interest points. The 'Interest Point 3: Scav wt pt' is selected. Below the dropdown, the 'Prescription' table is visible, showing the prescribed dose and number of fractions.

Rx ID	Rx Site	Prescribe To	Rx Dose (cGy)	Number of Fractions	Fractional Dose (cGy)
Physician's Intent	A	SCF	Interest Point 3: Scav wt pt	5000.0	25

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Remove Dosimetric Criterion

Remove All Dosimetric Criteria

- Minimum Dose
- Maximum Dose
- Mean Dose (Lower Limit)
- Mean Dose (Upper Limit)
- Minimum Dose Received by Relative Volume
- Minimum Dose Received by Absolute Volume
- Maximum Dose Received by Relative Volume
- Maximum Dose Received by Absolute Volume
- Minimum Relative Volume That Receives Dose
- Minimum Absolute Volume That Receives Dose
- Maximum Relative Volume That Receives Dose
- Maximum Absolute Volume That Receives Dose

Add Dosimetric Criterion

Structure: RECTUM

Dose (cGy): 6000.0

Maximum volume (%): 30.00

Volume tolerance (%): 0.00

Description:
The volume that receives 6000 cGy is smaller than 30 % with a tolerance of 0 %.

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

MonProstate
QA
MonProstate-QA1
Plan1
ProstVMAT


Load
Load into Planning
Load into Plan Review
Unload
New Bias Dose Plan
Delete Plans
CT1: Sum Plans
CT1: MonProstate - ProstVMAT
CT1: ProstVMAT - MonProstate

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
Structures@USS139QYMH2 - [002445TRN, Head_and_Neck^TG244, HN, VMAT3]

Name	Visi...	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlin...
BODY	<input checked="" type="checkbox"/>	12603.890	External	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
BODY-PTV56	<input checked="" type="checkbox"/>	13904.708	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
BRAIN	<input checked="" type="checkbox"/>	822.461	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
BRAINSTEM	<input checked="" type="checkbox"/>	31.464	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
CTV56	<input checked="" type="checkbox"/>	827.512	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
CTV63	<input checked="" type="checkbox"/>	351.814	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
CTV70	<input checked="" type="checkbox"/>	127.748	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
CTV56-CTV63	<input checked="" type="checkbox"/>	518.418	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
CTV63-CTV70	<input checked="" type="checkbox"/>	346.216	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
GTV70	<input checked="" type="checkbox"/>	27.599	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
Headrest	<input checked="" type="checkbox"/>	104.582	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>
LARYNX	<input checked="" type="checkbox"/>	14.262	Internal	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>



Structures Prescription Beams IMRT Constraints Dose Reference Points DVH Statistics

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




The screenshot shows the Elekta treatment planning system interface. The left sidebar contains a menu with 'Treatment Couch' highlighted. The main window displays a 'Save Treatment Couch As' dialog box with 'Couch Name: VersaHD' and 'Save' and 'Cancel' buttons. A table on the right lists the couch properties:

ble	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlin
23288.528		Couch	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.600	<input checked="" type="checkbox"/>
18717.946		Couch	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.050	<input checked="" type="checkbox"/>
259.139		Couch	<input type="checkbox"/>	<input type="checkbox"/>		<input checked="" type="checkbox"/>

Annotations on the screenshot include:

- ① Save as Treatment Couchをクリック (Click Save as Treatment Couch)
- ② 名前を入力してSave (コピーを作成) (Enter name and Save (create copy))
- ③ 作成したCouchをロード (Load the created couch)
- ④ Force EDにチェック、Relative EDに数値を入力しSave (Check Force ED, enter value for Relative ED and Save)

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

Structures					
View:	Contoured	All	Layers	Adapt Setup	
Name	Color	Force ED	Fill ED	Relative ED	
Foam Core	Green	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.050	
Carbon Fiber	Red	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.600	
Metal	Yellow	<input type="checkbox"/>	<input type="checkbox"/>		

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
Beams										General	Geometry	Treatment Aids	Setup Beams
Delete Parent Beams													
Beam	Description	Wedge ID	Angle	Orient	Port	MLC	Applicator ID	Bolus	SBD (cm)	Couch			
1	Left ear	No Wedge					6x6	5mmBol	99.01				

Bolus	SBD (cm)
Bolus	100.00

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Beams										General	Geometry	Treatment Aids	Setup Beams
Delete Parent Beams													
Beam	Description	Field ID	Visible	Delivery	Treatment ...	Modality	Algorithm	Energy	MU / Fx	Setup	SSD (cm)		
1	Left ear	1-1	<input checked="" type="checkbox"/>	3D	AMCInfin	Electron	Monte Carlo	6.0 MeV	440.18	SSD	100.99		

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Prescription													
Prescription													
Add Rx Delete Rx													
Rx ID Rx Site Prescribe To Rx Dose (cGy) Number of Fractions Fractional Dose (cGy)													
Physician's Intent A Lt Ear Depth of Beam 4 0.31 cm 5000.0 25 200.0													

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Beams										General	Geometry	Treatment Aids	Setup Beams
Delete Parent Beams													
Beam	Description	Field ID	Visible	Delivery	Treatment ...	Modality	Algorithm	Energy	MU / Fx	Setup	SSD (cm)	Bolus	SBD (cm)
1	Left ear	1-1	<input checked="" type="checkbox"/>	3D	AMCInfin	Electron	Monte Carlo	6.0 MeV	440.18	SSD	100.99	Bolus	100.00

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

Prescription Segments

Add Rx Delete Rx

Rx ID	Rx Site	Prescription
▼ Physician's Intent	A	Center of CTV SCV

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

Prescription Segments

Add Rx

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Rx ID	Rx Site	Prescription
▼ Physician's Intent	B	Center of CTV Chestwall

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
- ①Target →Target Penalty or Target EUD
- ②Patient→Quadratic Overdose (or Conformality)

IMRT Constraints											
↑ ↓ Pareto Constrained IMRT Parameters											
Structure	Cost Function	Delete	Enabled	Status	Manual	Weight	Reference Dose (cGy)	Multicriterial	Power Law Exponent	Shrink Margin (...)	Isoconstraint
PTV	Target Penalty		<input checked="" type="checkbox"/>	On	<input type="checkbox"/>	1.00					7000.0
patient	Quadratic Overdose		<input checked="" type="checkbox"/>	On	<input type="checkbox"/>	0.01	7000.0	<input type="checkbox"/>		0.00	100.0
<click to add a new structure>											
Isoeffect Relative Impact											
											6394.8
											1.1

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Equivalent Uniform Dose
= 等価均一線量

Target Penalty@DESKTOP-N4TCLM8 - (PROSTATE, Fusion^Prostate, CT...

Required Parameters

Prescription (cGy): 7000.0

Minimum Volume (%): 95.00

Optional Physical Parameters

Surface Margin: ☒

Close

Target Penalty

Target EUD@DESKTOP-N4TCLM8 - (PROSTATE, Fusion^Prostate, CT1, ...

Required Parameters

Prescription (cGy): 7000.0

Cell Sensitivity: 0.50

Optional Physical Parameters

Surface Margin: ☒

Close

Target EUD

Cost Function – Quadratic Overdose

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Quadratic Overdose@DESKTOP-N4TCLM8 - [PROSTATE, Fusion^ Prosta...

Required Parameters

Maximum Dose (cGy): 7000.0

RMS Dose Excess (cGy): 100.0


Optional Physical Parameters

Multicriterial: ☐

Shrink Structures

Structure Name	Include	Margin (cm)
PTV	<input checked="" type="checkbox"/>	0.00

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Properties

Remove Structure

Target EUD
Target Penalty
Quadratic Overdose
Parallel
Serial
Overdose DVH
Underdose DVH
Maximum Dose
Quadratic Underdose
Conformality

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IMRT Constraints

Pareto Constrained IMRT Parameters

Structure	Cost Function	Del
PTV	Target Penalty	
patient	Quadratic Overdose	

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Structures Prescription Beams IMRT Constraints Dose Reference Point

Target Penalty@DESKTOP-N4TCLM8 - [PROSTATE, Fusion^ Prostate, CT...

Required Parameters

Prescription (cGy): 7000.0

Minimum Volume (%): 95.00

Optional Physical Parameters

Surface Margin: ☒

Close

IMRT Constraintsの編集

Cost Functionの削除


Cost Function	Delete	Enabled
Target Penalty	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Quadratic Overdose	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Conformality	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

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

Grid Settings

Grid Spacing (cm): 0.30

Calculate Dose Deposition to: Medium

Store and display dose in couch structures: ☐

Grid Settings changes will be applied to ALL Rx IDs.

Algorithm Settings

Algorithm: Monte Carlo Photon

Statistical Uncertainty (%): 1.00

☐ Per Control Point ☒ Per Calculation

Sequencing Parameters: Dyn. Conformal Arc

☒ Constant Dose Rate

☐ Segment Shape Optimization

Selecting SSO will result in changes to the apertures defined by DCAT.

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Rx ID	Rx Site	Prescribe To	Rx Dose (cGy)	Number of Fractions	Fractional Dose (cGy)
▼ Physician's Intent	A	Plan Isocenter	7000.0	35	200.0

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Rescale 7000.0 cGy to...

Rx ID	Rx Site	Prescribe To	Rx Dose (cGy)	Number of Fractions	Fractional Dose (cGy)
▼ Physician's Intent	A	Plan Isocenter	7000.0	35	200.0

Actual Dose = 7000.0 cGy

Rescale 7000.0 cGy to dose at point Plan Isocenter

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

QA Clinic: MonacoPhantom: 30x30x30

Studyset Orientation(30x30x30): Head First Supine

Treatment Plan Orientation(Plan1): Head First Supine

Select Studyset Orientation for QA Plan: ☒ Head First ☐ Feet First

Reset Beams to Nominal Angles

☒ Gantry ☒ Collimator ☒ Couch

Calc Vol Grid Spacing (cm): 0.20

Algorithm: Collapsed Cone

Force entire volume to be treated as water: ☒

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
QA Planの作成

Interest Points & Markers@DESKTOP-N4TCLM8 - [MonacoPhantom, 30x30x30, Monac...

New Interest Point
New Marker
Jump to Point
Delete

ID	X (cm)	Y (cm)	Z (cm)	Description
I1	0.00	0.00	15.00	SSD
I2	0.00	0.00	10.00	5cm
I3	0.00	0.00	5.00	10cm
I4	0.00	0.00	0.00	ISO

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
Plan1-QA1


Attention: Unit not approved for Clinical use: Test Unit

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

The screenshot shows the 'Setup Options' window with the following settings:

- DICOM Settings...: Manage and edit DICOM import and export settings.
- CT-to-ED Files...: Manage and edit CT-to-ED conversion files.

The 'CT-to-ED Conversion Files' window shows a table of CT and ED values:

CT	ED
1000	0.001
-773	0.190
-516	0.490
-72	0.950
-34	0.980
-4	1.000
42	1.040
49	1.050
238	1.120
951	1.510

A graph on the right shows the relationship between CT and ED values, with a curve passing through the data points.

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The screenshot displays the Elekta software interface with the 'CT-to-ED Conversion Files' dialog box open. The 'Clinic' dropdown is set to '0~Clinic' and the 'Name' dropdown is set to 'DICOM3.Legacy120KV'. The 'Set As Default' checkbox is checked. A graph on the right shows the conversion curve for CT to ED values.

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Confidential and Proprietary Information, © 2021 Elekta, Inc. All rights reserved.