AQUA

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
The screenshot displays the AQUA ACUMYN interface. The top navigation bar includes a user profile (aquaphysicist) and menu items: Tests, Scheduling, Equipment, Compliance, and System Data. The main navigation bar features the AQUA logo, ACUMYN, and a 'To-Do List' button (highlighted in red) along with icons for Dashboard, Status, Library, Events, and Review.

The 'To-Do List' section contains filter options for 'Machine Type' (Linear Accelerator) and 'Machine' (Infinity TPI). Below these are tabs for 'Linear Accelerator' and 'Infinity TPI (Elekta Infinity)', with a status indicator 'Working - Clinical'.

The 'Test List' section includes a dropdown menu set to '-- All --' and a 'Show Last Run' checkbox. The table below lists various tests with their performance and compliance status.

| Test Name | Period | Performance | Compliance Status |
|---|--------|-------------|-------------------|
| KV-MV Iso Check | Daily | Green | Red DUE |
| Therapist Daily All Linacs | Daily | Green | Red DUE |
| BeamChecker Output Flatness Symmetry (18 MV, 0 Deg) | Daily | Green | Green |
| BeamChecker Output Flatness Symmetry (6 MV, 0 Deg) | Daily | Green | Red DUE |
| DQA3 Lite Output Flatness Symmetry (18 MV, 0 Deg) | Daily | Yellow | Red DUE |
| DQA3 Lite Output Flatness Symmetry (6 MV, 0 Deg) | Daily | Green | Red DUE |
| Flatness and Symmetry ICPProfiler (6 MV) | Daily | Green | Red DUE |

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The screenshot displays the AQUA interface for a Linear Accelerator. At the top, there are two dropdown menus: 'Machine Type' (set to '-- All --') and 'Machine' (set to 'VersaHD A'). Below these, the machine is identified as 'Linear Accelerator VersaHD A (Elekta VersaHD)' with a status of 'Working - Non-Clinical'. A red banner indicates 'MANDATORY TESTS INCOMPLETE'. The 'Test List' section shows a table of tests with their respective periods, performance, and compliance statuses.

| Test Name | Period | Performance | Compliance Status |
|--|----------|-------------|-------------------|
| Mandatory Safety CNSC - Elekta | Daily | Green | Red DUE |
| Daily Output Constancy (8 MV) | Daily | Grey | Red DUE |
| Field Coincidence Jaw and MLC (8 MV, 10 cm X, 10 cm Y) | Monthly | Red | Green |
| Flatness and Symmetry (8 MV) | Monthly | Green | Yellow 27 DAYS |
| Monthly Output Constancy Baseline (8 MV) | Monthly | Green | Yellow 27 DAYS |
| Output Constancy (8 MV) | Monthly | Yellow | Green |
| IAEA 398 Photons (8 MV) | Annually | Green | Grey 272 DAYS |
| Photon RDF (8 MV) | Annually | Green | Grey 272 DAYS |

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| Period | Performance | Compliance Status |
|----------|---|--|
| Daily |  |  DUE |
| Daily |  |  DUE |
| Monthly |  |  |
| Monthly |  |  28 DAYS |
| Monthly |  |  28 DAYS |
| Monthly |  |  |
| Annually |  |  273 DAYS |
| Annually |  |  273 DAYS |
| Annually |  |  273 DAYS |

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| Period | Performance | Compliance Status |
|---------|---|---|
| Daily |  |  DUE |
| Daily |  |  DUE |
| Daily |  |  DUE |
| Daily |  |  DUE |
| Monthly |  |  |
| Monthly |  |  11 DAYS |
| Monthly |  |  11 DAYS |
| Monthly |  |  11 DAYS |
| Monthly |  |  11 DAYS |

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
Binary test

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To-Do List

A screenshot of the AQUA To-Do List form. The form is light gray and contains two dropdown menus. The first dropdown menu is labeled 'Machine Type' and has 'Linear Accelerator' selected. The second dropdown menu is labeled 'Machine' and has 'VersaHD A' selected. Both dropdown menus are highlighted with a red box.

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Linear Accelerator VersaHD A (Elekta VersaHD) Working - Non-Clinical

Machine Compliance (CNSC Mandatory Tests)

MANDATORY TESTS INCOMPLETE

Test List

-- All -- Show Last Run

| Test Name | Period | Performance |
|--------------------------------------|--------|-------------|
| Light Field Size Test | Daily | Green |
| Therapist Daily All Linacs | Daily | Green |
| Therapists Daily Tests - Elekta All | Daily | Green |
| Gantry Angle Readout - Elekta | Weekly | Green |
| Therapist Collimator Readout Agility | Weekly | Green |

Test Name

- Couch Isocenter
- Isocenter Runout Gantry
- Isocenter Runout Gantry Short
- Isocenter Runout Gantry Short LinacConnect (8 MV)
- Mandatory Safety CNSC - Elekta**
- MLC Leaf and Jaw Position (8 MV)

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Test Form

Intercom PASS FAIL

Procedure:
1- Check that the intercom is on.
2- Verify that noise from inside the room, i.e. machine noise during beam on, a radio, etc. can be clearly heard through the intercom. Adjust the volume control if necessary.
Specifications:
The test fails if the intercom is not functioning

Video Monitor PASS FAIL

Procedure:
1- Switch on the TV remote control unit and the video monitor.
2- Check that all positioning and focus controls work.
3- Check that the picture quality is adequate to detect any patient motion.
Specifications:
The test fails if the video monitor system is not functioning

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Test Form

| | | | |
|---|---------------------------------------|---------------------------------------|----------------------------|
| <input type="checkbox"/> Intercom | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Video Monitor | <input type="radio"/> PASS | <input checked="" type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Wall Oxygen | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Wall Suction | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Couch Locks | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Motion Enable Switches | <input type="radio"/> PASS | <input checked="" type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> MLC Leaf Test Pattern | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Cross Hair Evaluation | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |
| <input type="checkbox"/> Vacloks | <input checked="" type="radio"/> PASS | <input type="radio"/> FAIL | <input type="radio"/> SKIP |

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


| Test Name | Period | Performance | Compliance Status |
|-------------------------------------|--------|---|---|
| KV-MV Iso Check | Daily |  |  DUE |
| Light Field Size Test | Daily |  |  DUE |
| Therapist Daily All Linacs | Daily |  |  |
| Therapists Daily Tests - Elekta All | Daily |  |  DUE |

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| Test Name | Period | Performance | Compliance Status |
|----------------------------|--------|---|---|
| KV-MV Iso Check | Daily |  |  DUE |
| Light Field Size Test | Daily |  |  DUE |
| Therapist Daily All Linacs | Daily |  |  |

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Machine Type: -- All -- Machine: VersaHD A

Linear Accelerator: VersaHD A (Elekta VersaHD) Working - Clinical

Machine Compliance (CNSC Mandatory Tests)

+ ALL MANDATORY TESTS COMPLETE

Test List

-- All -- Show Last Run

| Test Name | Period | Performance | Compliance Status |
|--|---------|-------------|-------------------|
| Mandatory Safety CNSC - Elekta | Daily | — | — |
| KV-MV Iso Check | Daily | — | — |
| Output ConstancyPMH (6 MV) | Daily | — | — |
| Field Coincidence Jaw and MLC (6 MV, 10 cm X, 10 cm Y) | Monthly | — | — |

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The screenshot displays two sections of the manual, each with a red box around its title:

- Procedure**
 - Image Acquisition**
 1. Put the Linac in Clinical mode.
 2. In Mosaic, open patient Monthly QA, \$\$##, where ##\$\$ is the machine ID, i.e. ES05. Load the 6 MV @fieldsize field and bring the Linac to a Beam Ready state.
 3. On the Iview computer
 - a. Log on to Iview.
 - i. Click the "other" radio button. Type in, or navigate to, "D:\Qadatabase" in the text box, as shown in Figure 1.
 - ii. Press "OK".
 4. Open up the patient you wish to acquire
 - i. Select the test you are performing from the Patient list in the bottom left hand corner of the screen.
 - ii. Select either "Jaw-defined" or "MLC-defined" from the treatment listing, depending on which symmetry test you wish to perform
 - iii. Highlight the field you wish to deliver
 5. Enable imager for acquisition
 - i. Click on the Single Exposure icon found under the field list box. A window will appear saying "Ready to acquire port image."
 6. Extend the PI such that the entire field is inside the active area of the imaging panel.
 - i. Rotate the couch around the pedestal so it will not interfere with the imager's travel or image acquisition
 7. Set-up the light field templates as shown in figure 2. Line up the BBs of the cross-hair template with the cross-hair shadow on the EPID. Also line up each of the light field corner etchings with the boundary of each light field corner.
 8. Deliver the treatment beam. The image should appear on the Iview screen.
 - Image Export**
 9. Immediately after you have take an image, click on the image Export icon
 10. An export window will appear.
 - i. Ensure the "This image only" radio button is selected.
 - ii. Click on the DICOM radio button.
 - iii. Select the "Witch_QA" directory from the pull-down menu.
 - iv. Click "OK".
 - Image Analysis**
 11. Log onto any computer.
 12. Open the eQA program, Start Menu->Citrix Applications->QA Tools->eQA.
 13. On the main menu screen of eQA, click on the Analyze QA image button, shown in Figure 3.
 14. Choose the image you wish to analyze.
 - i. The default folder is \\Witch\Dicom\QA. Click on the folder which contains the image you want to analyze (example U10_Symmetry_6x).
 - ii. The file names in this folder are standard Dicom export names from Iview One aspect of Iview's Dicom naming convention is the files appear in chronological order. Therefore the last image exported, will be the last image to appear in the list in the open dialog box.
 15. Select the image file to be analyzed and press the open button.
 16. The program will display a window asking you if you want to perform the "Symmetry test". Click on "Yes"
 17. Choose the beam energy used to acquire the image.
 18. Once the analysis is complete the screen will look similar to Figure 4.
 - i. Accept or Discard these results
 - a) If you click accept "æ" a summary of the results of the test will be saved, to be emailed to the staff responsible for maintaining the machine.
 - b) If you click discard "æ" the results are discarded, a summary of the tests will not be saved.
- Specification**
 1. The test pass if:
All four jaws are within 0.1 cm of the nominal jaw position.
 2. The test fails if:
If at least one jaw differs from the nominal jaw position by more than 0.1 cm

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Coincidence Baseline(cm): 0.0

Coincidence Tolerance(cm): 0.2

| Jaw | Optical Field (cm) | Radiation Field (cm) |
|-----|--------------------------------|--------------------------------|
| Y2 | <input type="text" value="5"/> | <input type="text" value="5"/> |
| Y1 | <input type="text" value="5"/> | <input type="text" value="5"/> |
| X2 | <input type="text" value="5"/> | <input type="text" value="5"/> |
| X1 | <input type="text" value="5"/> | <input type="text" value="5"/> |

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Coincidence Baseline(cm): 0.0

Coincidence Tolerance(cm): 0.2

| Jaw | Optical Field (cm) | Radiation Field (cm) | Difference (cm) |
|-----|--------------------|----------------------|-----------------|
| Y2 | 5 | 5 | 0.0 |
| Y1 | 5 | 5 | 0.0 |
| X2 | 5 | 5 | 0.0 |
| X1 | 5 | 5 | 0.0 |

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To-Do List

Machine Type: -- All --
Machine: Infinity IB

| Test Name | Period | Performance | Compliance Status |
|---|--------|-------------|-------------------|
| Couch Isocenter | Daily | | DUE |
| Isocenter Runout Gantry | Daily | | DUE |
| Isocenter Runout Gantry Short | Daily | | |
| Isocenter Runout Gantry Short XVI | Daily | | |
| Isocenter Runout Gantry XVI | Daily | | DUE |
| MLC Leaf and Jaw Position (6 MV) | Daily | | DUE |
| MLC Leaf and Jaw Position LinacConnect (6 MV) | Daily | | DUE |

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PROCEED SKIP BEAM DELIVERY

GO **CANCEL**

| | |
|---|----------|
| Active Contour for Open Field (Coll = 180°) | Complete |
| Extended Leaf Detection (Coll = 180°) | Complete |
| Leaf Position Detection (Coll = 180°) | Complete |
| Jaw Position Detection (Coll = 180°) | Complete |
| Isocenter Runout Collimator | PASS |
| Merge Leaf Position - Results | PASS |

Y1 MLC BANK **Y2 MLC BANK**

| Leaf# | -80 | -20 | 20 | 60 | 100 |
|-------|--------|--------|-------|-------|--------|
| 1 | -80.23 | -20.40 | 19.70 | 59.88 | 100.03 |
| 2 | -80.29 | -20.40 | 19.78 | 60.03 | 100.06 |
| 3 | -80.24 | -20.41 | 19.85 | 59.95 | 100.04 |
| 4 | -80.34 | -20.35 | 19.85 | 59.98 | 100.03 |
| 5 | -80.25 | -20.29 | 19.93 | 59.99 | 100.01 |
| 6 | -80.28 | -20.31 | 19.88 | 59.94 | 99.95 |
| 7 | -80.13 | -20.29 | 19.84 | 60.01 | 100.01 |

SAVE **CANCEL**

26 | Focus where it n

Restricted Information and Basic Personal Data

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Dashboard

| Machine | Effective Date | Machine Status | Mandatory | Compliance | Performance |
|------------|----------------|---|---|---|---|
| VersaHD A | Oct, 31 2016 |  Working - Non-Clinical |  |  |  |
| Infinity A | Oct, 31 2016 |  QA |  |  |  |
| Synergy B | Oct, 31 2016 |  No Change |  |  |  |
| Infinity C | Oct, 31 2016 |  QA |  |  |  |
| Synergy S | Oct, 31 2016 |  Working - Non-Clinical |  |  |  |
| Clinac IX | Oct, 31 2016 |  Out Of Service |  |  |  |

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The screenshot shows the AQUA dashboard interface. At the top, there is a navigation bar with the user name 'aquaphysicist' and menu items: Tests, Scheduling, Equipment, Compliance, System Data. Below this is the AQUA logo and 'AQUA TRAINING' text. A secondary navigation bar contains icons for 'To-Do List', 'Dashboard' (highlighted in red), 'Status', 'Library', 'Events', and 'Review'.

The main content area is titled 'Dashboard' and features a machine detail card. The card displays the following information:

| Name | Vendor | Model | Machine Status |
|-----------|--------|---------|------------------------|
| VersaHD A | Elekta | VersaHD | Working - Non-Clinical |

Below the card are two buttons: 'TEST VIEW' and 'BACK TO DASHBOARD'.

Below the card are two summary tables:

Mandatory

| | Total | Complete | Warn | Late | Pass | Fail |
|-----|-------|----------|------|------|------|------|
| All | 8 | 0 | 0 | 8 | 0 | 0 |


Compliance

| | Total | Complete | Warn | Late |
|---------|-------|----------|------|------|
| Daily | 18 | 0 | 0 | 18 |
| Weekly | 2 | 0 | 0 | 0 |
| Monthly | 8 | 0 | 0 | 8 |
| Other | 0 | 0 | 0 | 0 |

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| Mandatory | | Total | Complete | Warn | Late | Pass | Fail |
|---|-----|-------|----------|------|------|------|------|
|  | All | 8 | 0 | 0 | 8 | 0 | 0 |

| Compliance | | Total | Complete | Warn | Late |
|---|---------|-------|----------|------|------|
|  | Daily | 5 | 1 | 0 | 4 |
|  | Weekly | 1 | 0 | 0 | 0 |
|  | Monthly | 1 | 1 | 0 | 0 |
|  | Other | 0 | 0 | 0 | 0 |

| Performance | | Pass | Fail |
|---|---------|------|------|
|  | Daily | 1 | 0 |
|  | Weekly | 0 | 0 |
|  | Monthly | 1 | 0 |
|  | Other | 0 | 0 |

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| | Tests done? | Pass / Fail? | Indicator |
|------------------|-------------------|---------------------|-----------|
| Mandatory | yes | pass | green |
| | yes | (at least one) fail | red |
| | no | - | red |
| | No test scheduled | - | white |

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| Compliance | Tests done? | Indicator |
|------------|---|-----------|
| | yes | green |
| | No (because there were no tests scheduled) | white |
| | No (but it hasn't entered the near-late period yet) | white |
| | No (and it entered the near-late period) | yellow |
| | No (and it's due or period missed) | red |

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
| | Pass / Fail? | Indicator |
|--|---|-----------|
| Performance (of the last test run) | pass | green |
| | Pass (but exceeded the warning tolerance) | yellow |
| | fail | red |

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[SUMMARY VIEW](#) [BACK TO DASHBOARD](#)

Show Only Issues Hide Warn/Late Timestamp

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The screenshot displays the AQUA Status page. The top navigation bar includes the user 'aquaphysicist' and menu items: Tests, Scheduling, Equipment, Compliance, System Data. The main navigation bar features the AQUA logo, 'AQUA TRAINING', and icons for To-Do List, Dashboard, Status (highlighted in red), Library, Events, and Review.

Status

Set Status

Machine: Infinity TPI

Current Status: Working - Clinical

New Status: Working - Clinical

Comments:

UPDATE STATUS

Event List

| Timestamp | Type | Machine Status | Comments |
|--------------------|---------------------------|------------------------|----------------------|
| Oct, 31 2016 15:19 | Begin clinical treatments | Working - Clinical | |
| Oct, 24 2016 15:32 | Status Change | Working - Non-Clinical | QA measurement start |
| Oct, 24 2016 14:56 | Begin clinical treatments | Working - Clinical | Working Clinical |

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
Comments

New Status

- Select --
- No Change
- Working - Clinical
- Working - Non-Clinical
- Working - Disabled
- Out Of Service
- Scheduled Down Time
- QA
- Machine Off
- Research

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Set Status

Machine: VersaHD A

Current Status: Working - Non-Clinical

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
The screenshot shows the AQUA web interface. The top navigation bar includes the user name 'aquaphysicist' and menu items: Tests, Scheduling, Equipment, Compliance, System Data. The main navigation bar features the AQUA logo, 'AQUA TRAINING', and icons for To-Do List, Dashboard, Status, Library (highlighted in red), Events, and Review.

The 'Library' section contains search filters:

- Test Name: Input field with 'Output' entered.
- Machine: Dropdown menu with 'Infinity TPI' selected.
- Time Period: Dropdown menu with '(None)' selected.
- Radio buttons for 'Basic' and 'Scheduled' (selected).
- A blue 'SEARCH' button.

Below the filters is a 'Tests' table:

| Name | Version |
|--------------------------------------|---------|
| BeamChecker Output Flatness Symmetry | 2 |
| DQA3 Lite Output Flatness Symmetry | 2 |
| Output Constancy PTW Auto | 4 |
| Output Constancy PTW with Linac | 13 |
| QuickCheck Output Flatness Symmetry | 2 |


At the bottom left of the table, it says '1 - 5 of 5'.

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The screenshot shows the AQUA user interface. At the top, there is a dark navigation bar with the user name 'aquamorningga' and a dropdown arrow. Below this is the AQUA logo and 'INSTITUTION NAME'. To the right of the logo are several navigation buttons: 'To-Do List', 'Dashboard', 'Status', 'Library' (highlighted in red), 'Events', and 'Review'. Below the navigation bar, the word 'Library' is displayed. The main content area features a search form with a text input field labeled 'Test Name', which is highlighted with a red border. Below the input field is a blue 'SEARCH' button. To the right of the input field are two radio buttons: 'Basic' (selected) and 'Scheduled'.

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
| Tests | |
|-------------------------------------|--|
| Name | |
| Daily Calibration/Flatness | |
| Daily Output Constancy | |
| Daily Safety Tests MRL | |
| Therapist Daily All Linacs | |
| Therapists Daily Tests - Elekta All | |
| Therapists Daily Tests - Varian All | |

Showing 1 - 6 of 6

| Test Instances | |
|---|--|
| Name | |
| Therapist Daily All Linacs - Clinac IX | |
| Therapist Daily All Linacs - Infinity A | |
| Therapist Daily All Linacs - Infinity B | |
| Therapist Daily All Linacs - Infinity TPI | |
| Therapist Daily All Linacs - Synergy S | |
| Therapist Daily All Linacs - VersaHD A | |

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
Test Form

| | |
|---|---|
| <input type="checkbox"/> Intercom | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Video Monitor | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Wall Oxygen | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Wall Suction | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Couch Locks | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Motion Enable Switches | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> MLC Leaf Test Pattern | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Cross Hair Evaluation | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |
| <input type="checkbox"/> Vacloks | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> SKIP |

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The screenshot displays the AQUA web interface. At the top, a navigation bar includes the user name 'aquaphysicist', menu items like 'Tests', 'Scheduling', 'Equipment', 'Compliance', and 'System Data', and the AQUA logo with 'AQUA TRAINING'. A secondary navigation bar contains icons for 'To-Do List', 'Dashboard', 'Status', 'Library', 'Events' (highlighted in red), and 'Review'.

Below the navigation, the 'Events' section is active, with the 'Event List' tab selected and highlighted with a red box. The 'Clinical Dates' tab is also visible.

The main content area features a search filter with three dropdown menus: 'Machine' (set to 'VersaHD A'), 'Event Type' (set to '-- All --'), and 'Date Range' (set to 'Oct, 01 2015' to 'Oct, 31 2016'). A blue 'SEARCH' button and a grey 'NEW EVENT' button are positioned below the filters.

Below the search area is a table with the following data:


| Timestamp | Type | Machine Status | Comments |
|--------------------|---------------|------------------------|-----------------|
| Oct, 24 2016 15:37 | Annual QC | Working - Non-Clinical | Start annual QC |
| Jul, 25 2016 07:46 | Status Change | Working - Clinical | |
| Jul, 25 2016 07:28 | Status Change | Working - Clinical | |
| Jul, 25 2016 06:32 | Status Change | Working - Non-Clinical | |

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The screenshot shows the AQUA interface with the following elements:

- Header: AQUA ACUMYN INC. with navigation icons for To-Do List, Dashboard, Status, Library, Events (highlighted), and Review.
- Sub-headers: Event List and Clinical Dates.
- Search filters: Machine (VersaHD A), Event Type (-- All --), and Date Range (02/28/2019 to 03/14/2019).
- Buttons: SEARCH and NEW EVENT.
- Table of event records:

| Timestamp | Type | Machine Status | Comments |
|------------------|---------------------------|------------------------|----------|
| 03/14/2019 09:42 | Begin clinical treatments | Working - Clinical | |
| 03/14/2019 09:38 | Status Change | Working - Non-Clinical | |
| 03/14/2019 09:34 | Begin clinical treatments | Working - Clinical | |
| 03/14/2019 08:36 | Begin QA | QA | |

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
The screenshot illustrates the process of creating a new event in the AQUA system. On the left, the 'NEW EVENT' button is highlighted with a red box. A large pink arrow points to the right, where the 'Event List > New Event' form is shown. The form fields for 'Machine' (Infinity Test1), 'Event Type' (Begin QA), and 'Machine Status' (QA) are highlighted with red boxes. The 'Comments' field contains 'Morning QA' and is also highlighted with a red box. The 'SAVE' button is highlighted with a red box.

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The screenshot displays the AQUA interface for searching clinical dates. The 'Machine' dropdown menu is set to 'Flexitron A', and the 'Date Range' is set to '04/02/2020 - 04/16/2020'. A blue 'SEARCH' button is highlighted. Below the search form, the 'Clinical Date List' table shows two entries for 'Flexitron A' with dates '04/16/2020' and '04/09/2020'. The table is summarized as follows:

| Name | Date |
|-------------|------------|
| Flexitron A | 04/16/2020 |
| Flexitron A | 04/09/2020 |

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

The screenshot shows a search interface with a 'Machine' dropdown menu set to 'Flexitron A' and a 'Date Range' of '04/02/2020 - 04/16/2020'. A blue 'SEARCH' button is visible. Below it, a grey 'NEW DATE' button is highlighted with a red box.


The screenshot shows the 'Clinical Date List > New Clinical Date' form. The 'Machine' dropdown is set to 'Flexitron A'. The 'Date' field contains '04/02/2020'. Below the form, a blue 'SAVE' button, a red 'DELETE' button, and a dark grey 'CANCEL' button are visible. The 'SAVE' button and the 'Date' field are highlighted with red boxes.

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Review

Reports Plot Data **Test Results**

Test
Gantry Angle Readout - Elekta Include Draft Tests

Test Instance
Gantry Angle Readout - Elekta - Infinity A

Date Range
05/01/2018 - 09/07/2018

Test Runs
 Non-Draft Draft

SEARCH No Plot Available

Test Run List

| Run ID | Date | Complete | Status | Tester | Revision Of Run |
|--------|------------|------------------|---|---------------|-----------------|
| 447572 | 09/03/2018 | 09/03/2018 07:44 | ■ PASS | aquaPhysicist | |
| 446516 | 08/30/2018 | 08/30/2018 07:37 | ■ PASS | aquaPhysicist | |
| 442794 | 08/23/2018 | 08/23/2018 07:35 | ■ PASS | aquaPhysicist | |
| 436453 | 08/12/2018 | 08/12/2018 07:48 | ■ PASS | aquaPhysicist | |

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Test Run List

| Run ID | Date | Complete | Status | Tester | Revision Of Run |
|--------|------------|------------------|--------|---------------|-----------------|
| 447572 | 09/03/2018 | 09/03/2018 07:44 | PASS | aquaPhysicist | |
| 446516 | 08/30/2018 | 08/30/2018 07:37 | PASS | aquaPhysicist | |
| 442794 | 08/23/2018 | 08/23/2018 07:35 | PASS | aquaPhysicist | |
| 436453 | 08/12/2018 | 08/12/2018 07:48 | PASS | aquaPhysicist | |

Review Test

Test Instance Gantry Angle Readout - Elekta - Infinity A | Machine Infinity A | Status PASS

Date Jul, 19 2016 | Complete Jul, 19 2016 07:46 | Test Name aquaPhysicist

[SHOW INPUT FORM](#) [BACK](#)

Test Results

+ Procedure

+ Specifications

Please enter the digital and mechanical readouts below.

| Angle | Digital | Digital Result |
|----------------|---------|----------------|
| 0 | 0 | PASS |
| 90 | 89.8 | PASS |
| 180 | 180.2 | PASS |
| 270 | 269.7 | PASS |
| Overall Status | | PASS |

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The image shows a two-part screenshot of the AQUA interface. The left part shows a test summary with a 'SHOW INPUT FORM' button highlighted in a red box. A large pink arrow points from this button to the right part of the screenshot. The right part shows the 'Input Form' page, where the 'SHOW INPUT FORM' button is now 'SHOW RESULT FORM'. Below this, there is a table for entering digital and mechanical readouts, with the table area highlighted in a red box.

Test Instance: Gantry Angle Readout - Elekta - Infinity A | Machine: Infinity A | Status: PASS

Date: 09/03/2018 | Complete: 09/03/2018 07:44 | Tester: aquaPhysicist

Procedure

Specification

SHOW RESULT FORM | BACK | REVISE TEST RUN

Input Form

Procedure

Specifications

Please enter the digital and mechanical readouts below.

| Angle | Digital |
|-------|---------|
| 0 | 0.1 |
| 90 | 90 |
| 180 | 179.5 |
| 270 | 269.5 |

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The screenshots illustrate the workflow for revising a test run:

- Initial State:** The 'REVISE TEST RUN' button is highlighted in red.
- Calculation:** The 'CALCULATE' button is highlighted in red. The digital readouts table is as follows:


| Angle | Digital |
|-------|---------|
| 0 | 1 |
| 90 | 89 |
| 180 | 182 |
| 270 | 271 |
- Final State:** The 'SAVE' button is highlighted in red. The digital readouts table now includes a 'Digital Result' column:

| Angle | Digital | Digital Result |
|-------|---------|----------------|
| 0 | 1 | PASS |
| 90 | 89 | PASS |
| 180 | 181 | PASS |
| 270 | 271 | PASS |

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
Include Draft Tests

Test Plot
[Empty dropdown]

Event Types
(None)
Return to clinical service
Planned servicing
Annual QC
Out of Service

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Include Draft Tests

Test Plot
Result

Event Types
Return to clinical service
Planned servicing
Annual QC
Out of Service
Output adjustment

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
Daily Safety Tests

| Test Name | Expected(#) | Passed(#) | Failed(#) | Compliance(%) |
|---|-------------|-----------|-----------|---------------|
| Audible LPO Alarm | 86 | 85 | 0 | 98.8 |
| Beam On Console Warning Light - KV | 86 | 73 | 0 | 84.9 |
| Beam On Treatment Room Warning Light - KV | 86 | 73 | 0 | 84.9 |
| Control Warning Lights | 86 | 85 | 0 | 98.8 |
| Couch Locks | 86 | 74 | 0 | 86 |
| EPID and OBI Touchguard | 86 | 73 | 0 | 84.9 |
| Intercom | 86 | 85 | 0 | 98.8 |
| KV Door Interlock | 86 | 73 | 0 | 84.9 |
| Last Person Out Button Time-Out | 86 | 85 | 0 | 98.8 |
| Motion Enable Switches | 86 | 74 | 0 | 86 |
| Radiation Monitor | 86 | 85 | 0 | 98.8 |
| Room Warning Lights | 86 | 85 | 0 | 98.8 |
| Run Beam | 86 | 73 | 0 | 84.9 |
| Treatment Door Panic Panel and Crashbar | 86 | 85 | 0 | 98.8 |
| Treatment Room Door Interlock | 86 | 85 | 0 | 98.8 |
| Treatment Stop Counter Varian | 86 | 85 | 0 | 98.8 |
| Video Monitor | 86 | 85 | 0 | 98.8 |
| Wall Oxygen | 86 | 74 | 0 | 86 |
| Wall Suction | 86 | 74 | 0 | 86 |

Daily QC Tests

| | | | | |
|------------------------------|----|----|---|------|
| Cross Hair Evaluation | 86 | 74 | 0 | 86 |
| Lasers After KV/MV Alignment | 86 | 72 | 0 | 83.7 |
| Light Field Size Test | 86 | 74 | 0 | 86 |
| MLC Leaf Test Pattern | 86 | 74 | 0 | 86 |
| ODI Readout | 86 | 73 | 0 | 84.9 |
| On-Board Imaging QA | 86 | 74 | 0 | 86 |

Weekly and Monthly QC Tests

| | | | | |
|-------------------------------|----|----|---|------|
| Flatness and Symmetry (8 MV) | 4 | 4 | 0 | 100 |
| Gantry Angle Readout - Varian | 17 | 15 | 0 | 88.2 |

Compliance Level: Requires Follow Up

Report Generated by: AQUA Demo

Thank you

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COPYRIGHT 2007-2017 ACUMYN
Acumyn Inc. - Clinac IX - SN: 1011