# AQUA 2.0 to 3.0 - Apps Upgrade Training

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# 2- Automated Image-based Daily QA test

## Output Uniformity and MLC test

### Test Results

Waiting for Image Export	Complete
Image Import	Complete
Field Center and Edge Detection (1)	Complete
Output (1)	PASS
Field Center and Edge Detection (2)	Complete
Output (2)	PASS
Field Center and Edge Detection (3)	Complete
Output (3)	PASS
Average Measurements	Complete
Mean Isocenter Detection	Complete
Active Contour for Open Field, Collimator = 0°	Complete
Extended Leaf Detection, Collimator = 0°	Complete
Leaf Position Detection, Collimator = 0°	Complete
Jaw Position Detection, Collimator = 0°	Complete
Isocenter Runout Collimator	PASS
Leaf Position - Results	PASS
Jaw Position - Results	PASS

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


# 6- Compliance streamlined

## Mandatory/Optional compliance added to the Scheduling page

The screenshot shows the 'Scheduling' page with a search bar and a 'New Schedule' link. The form includes dropdown menus for 'Test' (Output Uniformity), 'Test Instance' (Output Uniformity - VersaHD A (8 MV)), 'User or Group' (aquaPhysicist), and 'Time Period' (Daily). There is an 'Auto Schedule' checkbox. A highlighted pink box contains the 'Mandatory' and 'Optional' checkboxes, with 'Mandatory' selected. Below it, the 'Include in Compliance Report (Mandatory Section)' checkbox is also selected. 'SAVE' and 'CANCEL' buttons are at the bottom.

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The screenshot shows a dashboard with a table of test results. The table is divided into two sections: Mandatory and Compliance. The Mandatory section has columns for Total, Complete, Warn, Late, Pass, and Fail. The Compliance section has columns for Total, Complete, Warn, and Late. The table shows one mandatory test (All) and four compliance tests (Daily, Weekly, Monthly, Other).

Mandatory		Total	Complete	Warn	Late	Pass	Fail
	All	1	0	0	1	0	0

Compliance		Total	Complete	Warn	Late
	Daily	1	0	0	1
	Weekly	0	0	0	0
	Monthly	1	1	0	0
	Other	0	0	0	0

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
The screenshot shows the 'Review' page in the AQUA ELEKTA system. The 'Test Results' tab is selected. The search filters are set to 'Audiovisual Monitors Intercom' for the Test and 'Audiovisual Monitors Intercom - VersaHD A' for the Test Instance. The date range is from 12/01/2019 to 12/31/2019. The 'Test Runs' section shows 'Non-Draft' selected. The 'Test Run List' table has the following data:

Run ID	Date	Complete	Status	Tester	Revision Of Run	Hidden
60745	12/19/2019	12/19/2019 11:04	PASS	aquaPhysicist		<input checked="" type="checkbox"/>
60743	12/18/2019	12/18/2019 13:13	FAIL	aquaPhysicist		<input type="checkbox"/>
60737	12/18/2019	12/18/2019 11:13	PASS	aquaPhysicist		<input type="checkbox"/>

The 'Hidden' checkbox for the first row is highlighted with a red circle. A 'SAVE' button is visible at the bottom right of the table.

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The screenshot shows the AQUA ELEKTA Review interface. The top navigation bar includes 'To-Do List', 'Dashboard', 'Status', 'Library', 'Events', and 'Review'. The 'Review' section has tabs for 'Reports', 'Plot Data', 'Test Results', and 'Approvals'. Below the tabs are filters for 'Machine' (VersaHD A), 'Date Range' (12/12/2019 to 12/19/2019), 'Time Period' (-- All --), 'Test' (-- Select --), and 'Approval Status' (Not Approved, All). A 'SEARCH' button is present. The 'Test Run List' table has columns: Run ID, Test Instance Name, Complete, Status, Approval Time, Approver, and Approved. The 'Approved' column contains checkboxes, with the first one checked. 'SELECT ALL & SAVE' and 'SAVE' buttons are at the bottom right. The footer shows 'Showing 1 - 3 of 3'.

Run ID	Test Instance Name	Complete	Status	Approval Time	Approver	Approved
60745	Audiovisual Monitors Intercom	12/19/2019 11:04	PASS	12/19/2019 11:52	aquaPhysicist	<input checked="" type="checkbox"/>
60743	Audiovisual Monitors Intercom	12/18/2019 13:13	FAIL			<input type="checkbox"/>
60742	Beam On Indicator Room Warning Lights	12/18/2019 11:37	PASS			<input type="checkbox"/>

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## Equipment

- Machines
- Machine Types
- Machine Type Groups
- Machine Groups
- Machine Components

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

# Appendix A

## Conversion Factor Setup

The screenshot displays the 'Edit Test Instance' page for 'Output Uniformity - VersaHD A (8 MV)'. The 'Parameters' section is highlighted, showing a table with the following data:

Label	Name	Type	Value
Conversion factor for output calculation	conversion_factor	test	9.91289E-5
Delivery Energy	energy	equipment	8 MV
Fail Tolerance	failTolerance	test	3.0 %
Image Directory and MRN	directory	test	Output Uniformity
Invert Image	invert_image	test	1.0
Output Baseline	outputBaseline	test	20.0 A.U.
Output Fail Tolerance	outputFailTolerance	test	3.0 %
Output Warn Tolerance	outputWarnTolerance	test	2.0 %
Radial Beam Shape Baseline	radialFlatnessBaseline	test	3.0 %
Radial Symmetry Baseline	radialSymmetryBaseline	test	0.0 %
ROI Percentage	openFieldPct	test	80.0
Transverse Beam Shape Baseline	transverseFlatnessBaseline	test	3.0 %
Transverse Symmetry Baseline	transverseSymmetryBaseline	test	0.0 %
Warn Tolerance	warnTolerance	test	2.0 %

An inset window titled 'Parameters' provides a magnified view of the first row, with the value '9.91289E-5' highlighted in a pink box. An arrow points from this inset to the corresponding row in the main table.

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
The screenshot shows a 'Test Results' window with a table of 'Average Measurements'. The table has two columns: 'Item' and 'Average Measurement'. The first row is highlighted in pink and contains 'Output in 2x2 CM Square Region (A.U. = Arbitrary Unit)' with a value of '100.272 A.U.'. The second row is 'Output STD of the 2x2 CM Square Region (A.U. = Arbitrary Unit)' with a value of '0.022 A.U.'. A dashed arrow points from this table to the table above.

Output value from the average measurements:

Item	Average Measurement
Output in 2x2 CM Square Region (A.U. = Arbitrary Unit)	100.272 A.U.
Output STD of the 2x2 CM Square Region (A.U. = Arbitrary Unit)	0.022 A.U.

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