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
The screenshot displays the Monaco Scripting API Index. The left sidebar shows a tree view of the API structure, with 'ExitMonaco Method' selected. The main content area shows the details for 'MonacoApplication.ExitMonaco Method'. The title is 'MonacoApplication.ExitMonaco Method' with the Elekta logo. Below the title, there is a description: 'Exit Monaco. Once the system execute the method, it will pause 1 second to check whehter the "Save Changes" dialog box occurs. If the "Save Changes" dialog box does not occur, Monao exits; If the "Save Changes" dialog pops out, you can use the handle it with the parameters below. Default is SaveDlgOptions.DoNothing. if no parameter is used.' The namespace is 'Elekta.MonacoScripting.API' and the assembly is 'MonacoScripting (in MonacoScripting.dll) Version: 6.1.2.0'. Under the 'Syntax' section, the C# code is shown:

```
C#  
public bool ExitMonaco(  
    MonacoApplication. SaveDlgOptions option = MonacoApplication.SaveDlgOptions.DoNothing  
)
```

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
```
/* Add your script steps. For example,  
* 1. Launch and login Monaco  
* 2. load a patient  
* ....  
* Close the Patient  
* Exit Monaco */
```

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The screenshot shows a software interface with a left-hand navigation pane and a main content area. The navigation pane lists various methods under the 'MonacoApplication Class' and 'MonacoApplication Methods' categories. The main content area displays the details for the 'MonacoApplication.GetCurrentPatient Method'. It includes a description: 'Get current patient demographic.', the namespace 'Elekta.MonacoScripting.API', and the assembly 'MonacoScripting (in MonacoScripting.dll) Version: 6.1.2.0'. Under the 'Syntax' section, a C# code snippet is shown: `public PatientDemographic GetCurrentPatient();`. Below the syntax, the 'Return Value' is specified as 'Type: PatientDemographic' and 'The PatientDemographic object, or null if it fails'. A 'See Also' section is also present.

The screenshot shows a web-based API documentation page. On the left is a navigation tree with the following items: Copyright and Trademarks, Elekta MonacoScripting API, MonacoApplication Class, MonacoApplication Properties, MonacoApplication Methods (expanded), AdaptAnatomy Method, ClosePatient Method, DVHExport Method, ExecuteBatchOptimization Method, ExecuteCalculationForActiveRX Method, ExecuteCalculationForAIRX Method, ExecuteOptimizationPhase1 Method, ExecuteOptimizationPhase2 Method, ExecuteOptimizeWithCurrentSpotWeight Method, ExecuteOptimizeWithInitialSpotWeight Method, ExitMonaco Method, ExportAllCustomizedReports Method, ExportCustomizedReports Method, ExportDosePlane Method, ExportReport Method, FindPatient Method, **GetActivePlan Method** (highlighted), GetAssignCTToEDSettings Method, GetAutoMarginCreator Method, GetBeamsSpreadsheet Method, GetBeamVisibilitySettings Method, GetCalculationPropertiesSettings Method.

MonacoApplication.GetActivePlan Method

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