# Unity Kickoff

## Part 1

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Ring Out

# Personalized adaptive radiotherapy

## OFFLINE

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Patient Preparation

Simulation

Reference plan

## ONLINE

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Patient set up

Daily MR

Adapt to position

Rigid registration

Adapt to shape

Deformable registration

Fluence re-optimization or Segment re-optimization (WSO)

Motion monitoring

Aperture shape adaptation

Beam delivery

# Patient Preparation

# MOSAIQ Patient Registration

**MOSAIQ**

- ✓ Patient Name
- ✓ Patient ID
- ✓ Birth Date
- ✓ Gender

# Patient Preparation

**MOSAIQ**

# MR Initial Assessment

Assessments/MR Initial Assessment - MR#: 4DCIRS1 One, CIRS

Diagnosis: Malignant neoplasm other parts pancreas [C25.7] Stage: Course: No Courses  
Histology: MD: Crumble, Beverly

Flowsheet | Clinician Worksheet | Laboratory | Vital Signs | Assessments | Graphs | View: MR Initial Assessment

Date	8/29/17
Time	9:25 AM
Pacemaker	2
Implants, wires or other	2
Voice prosthesis	2
Last 6 mnths endoscope	2
Dentures/ dental plate/ br	2
Hearing aid	2
Metal particles in your ey	2
Shrapnel or bullets in bod	2
Breast expander w/ magn	2
Hydrocephalus/insulin/me	2
Surgery-head,brain,eyes,	3
Surgery in past 2 months	3
Tattoos or skin patches	3
Colored contact lenses	3
Pregnant	2
Weight (lb) lb	155
Comment	

Buttons: Close, Add, Change, Delete, Copy, Status, Course, MAR Sum

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Abdomen-HFS-Motion[btFFE RealTime\_12min.]

Abdomen-HFS-Pre[T2 3D Tra 3min]

Activity Association to Prescription

Activity	Prescription
Abdomen-HFS-Pre[T2 3D Tra 3min]	Pancreas 7,400.00 cGy / 37# 200.00 cGy

Image Acquisition

Modality: MRI Exam Card: Abdomen Acquisition: T2 3D Tra 3min Purpose: Plan Adaptation

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

# Simulation

## Record setup and reference mark in Site Setup Definition

MOSAIQ

Site Setup Definition - MR#: MRLProstate Prostate, Unity

Rx Site: training      Dose: ??????/7,800 cGy      Fractions: ??/39      Approved:

Site Setup

Patient Orientation: Head In, Supine      Machine: ElektaMRL      Created:      Approved:      Last Modified:

Setup

Coils and Patient References

Reference Mark	Index
25.0	

[Add Coil]

Accessories

Accessory	Index
Handgrip M	
Headrest M	14.0
Kneestep+block M	29

[Add device]

Photos/Diagrams

Setup Instructions

REMOVE HEARING AID DAILY

Scan Reference Offset is 1.0 cm away from 29.5, towards foot of bed

Couch

Site Setup Definition Will Be Added

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

# Personalized adaptive radiotherapy

## OFFLINE

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Patient Preparation


Simulation

Reference plan

## ONLINE

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Patient set up

Daily MR

Adapt to position

Rigid registration

Adapt to shape

Deformable registration

Fluence re-optimization or Segment re-optimization (WSO)

Motion monitoring

Aperture shape adaptation

Beam delivery

# Patient Set-up

Patient specific set-up details presented at console and in treatment room

MOSAIQ

Head In, Supine

Set: 20.5 cm  
Actual: 17.3 cm  
Table Not Ready

Patient Positioning Device & Reference Mark	Index
Headrest - Supine	
Wing Step	3
Knee Step	35
Feet Step	45
Reference Mark	12

Coil direction	Index	Direction
Coil frame	16	→ Bone
Anterior Coil		

Setup Photos 1/4

Check the checkbox to select two photos to be displayed in treatment room

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# Daily MR Checks

Check the form to ensure the patient is eligible for a daily treatment session

**MOSAIQ**

MR Daily Assessment - MR#: 12345 ATLANTIC, John

Record Date and Time: 10/11/2016 09:54:09

Has anything changed?

Pacemaker:

Implants, wires or other:

Metal particles in your eyes:

Shrapnel or bullets in body:

Breast expander w/ magnet:

Hydrocephalus/insulin/med pump:

New Surgery since last visit:

Tattoos or skin patches:

Colored contact lenses:

Pregnant:

Weight (kg):

Comment:

MR Daily Assessment Will Be Added

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Elekta

Unity

PHILIPS

# Personalized adaptive radiotherapy

## OFFLINE

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Patient Preparation

Simulation

Reference plan

## ONLINE

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Patient set up

Daily MR

Adapt to position

Rigid registration

Adapt to shape

Deformable registration

Fluence re-optimization or Segment re-optimization (WSO)


Motion monitoring

Aperture shape adaptation

Beam delivery


# Image Acquisition

**MOSAIQ**

Start imaging based on presets defined in Plan of Care

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# Personalized adaptive radiotherapy

## OFFLINE

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Patient Preparation

Simulation

Reference plan

## ONLINE

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Patient set up

Daily MR

Rigid registration


Adapt to shape

Adapt to position

Deformable registration

Fluence re-optimization or Segment re-optimization (WSO)

Aperture shape adaptation

Motion monitoring

Beam delivery


# Image Registration

Monaco

## Automatic initial registration of the anatomy of the day

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



# Adapt to Position

Choose between adapt to position or adapt to shape based on clinical need.

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
# Personalized adaptive radiotherapy

## OFFLINE

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Patient Preparation

Simulation

Reference plan

## ONLINE

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Patient set up

Daily MR

Rigid registration

Adapt to shape

Adapt to position

Deformable registration

Fluence re-optimization or Segment re-optimization (WSO)

Motion monitoring

Aperture shape adaptation

Beam delivery

# Delivery

## New plan is delivered and recorded with simultaneous motion monitoring

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Day 1 Classroom / Dept. Training	Day 2 Classroom / Dept. Training	Day 3 Bunker Training	Day 4 Bunker Training	Day 5 Bunker Training
<p><b>AM</b></p> <ul style="list-style-type: none"> <li>Introduction</li> <li>Safety Overview</li> <li>Geography Overview</li> <li>System Demo and Overview</li> </ul> <p><b>PM</b></p> <ul style="list-style-type: none"> <li>Patient Prep</li> <li>Clinical Administration</li> <li>Simulation</li> <li>Contouring Review</li> </ul>	<p><b>AM</b></p> <ul style="list-style-type: none"> <li>Initial Treatment Plan Creation</li> <li>Hands on Initial Planning</li> <li>Pre-Treatment Chart Prep</li> </ul> <p><b>PM</b></p> <ul style="list-style-type: none"> <li>Optimization Parameters</li> <li>Online Planning Options and Practical Adaptive Planning Sessions</li> </ul>	<p><b>AM</b></p> <ul style="list-style-type: none"> <li>Adapt-to-position Workflow with Motion Monitoring and Treatment Delivery</li> <li>Adapt-to-position Workflow with Hands-on Session</li> </ul> <p><b>PM</b></p> <ul style="list-style-type: none"> <li>Adapt-to-position Workflow with Hands-on Session</li> <li>Adapt-to-position Workflow with Verification Image</li> </ul>	<p><b>AM</b></p> <ul style="list-style-type: none"> <li>Adapt-to-shape Workflow Hands-on Session</li> </ul> <p><b>PM</b></p> <ul style="list-style-type: none"> <li>Adapt-to-shape Workflow Hands-on Session</li> <li>Treatment Session Interruptions Workflow Hands-on Session</li> </ul>	<p><b>AM</b></p> <ul style="list-style-type: none"> <li>Treatment Session Interruptions Workflow Hands-on Session</li> </ul> <p><b>PM</b></p> <ul style="list-style-type: none"> <li>Navigator Imaging Workflow</li> <li>Navigator Imaging Workflow Hands-on Session</li> <li>Imaging-only Workflow</li> <li>Imaging-only Hands-on Session</li> </ul>

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