Unity Physics Training

Unity Software System Overview

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

Unity Room Layout - Control Room

Software Components

Review of Unity Applications Software

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
The screenshot displays the Monaco software interface for treatment planning. The main window shows a cross-sectional view of a patient's head and neck with a color-coded dose distribution. A Monaco logo is visible in the top right corner. Below the main view, there is a table with treatment session details and a Monaco logo in the bottom right corner.

Tx	Field	Status	Meterset	Dose	E/M	Pattern	Site Setup
1			70.6 MU	24 cGy	7 X FFF StepNShoot 18 Control Poi		
2			45.2 MU	24 cGy	7 X FFF StepNShoot 12 Control Poi		
3			112.5 MU	3 cGy	7 X FFF StepNShoot 28 Control Poi		
4			73.3 MU	24 cGy	7 X FFF StepNShoot 18 Control Poi		
5			48.2 MU	20 cGy	7 X FFF StepNShoot 12 Control Poi		
6			58.6 MU	4 cGy	7 X FFF StepNShoot 12 Control Poi		
7			82.8 MU	3 cGy	7 X FFF StepNShoot 22 Control Poi		

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

Introduction


Computer Systems

Technical Room Computer Cabinet (TRCC)

Overview

- Clinical Administration
- Plan of Care
- MOSAIQ
- Monaco
- TSM

User Manual

Elekta Unity
Instructions for Use - Clinical Mode


Applications
Software
Components

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
The screenshot displays the 'Image Library > Acquisition > MRI' section of the software. It features a search bar with 'MRI' entered and a 'Quick Search' button. Below is a table with two columns: 'Exam Card' and 'Acquisition Name'. The table lists various MRI sequences and their corresponding acquisition names.

Exam Card	Acquisition Name
MRL Abdomen FF	SURVEY
MRL Abdomen FF	SURVEY_optional (OPT)
MRL Abdomen FF	T2 3D Tra
MRL Abdomen FF	T1 3D Tra
MRL Abdomen FF	T2 3D Tra Navigated (NAV)
MRL Abdomen FF	btFFE Cor RealTime
MRL Abdomen FF	btFFE Sag RealTime
MRL Abdomen FF	btFFE Cor Sag RealTime
MRL Abdomen FF	btFFE Sag Tra RealTime
MRL Abdomen FF	btFFE Cor Sag Tra RealTime
MRL Abdomen HF	SURVEY
MRL Abdomen HF	SURVEY_optional (OPT)
MRL Abdomen HF	T2 3D Tra
MRL Abdomen HF	T1 3D Tra
MRL Abdomen HF	T2 3D Tra Navigated (NAV)
MRL Abdomen HF	btFFE Cor RealTime

Unity Applications Software

Clinical Administration: Applying Activity Templates during the Online Workflow

The screenshot displays the Unity Applications Software interface for a treatment session. The patient information panel shows: Rectum, One; MR#: rectumb1; 12/12/1970; 47 yr; 154.1* Rectum, NOS Affirmed. The main workspace contains four image views (axial, sagittal, coronal, and another axial view) showing the rectum and surrounding structures. The Structure Visibility list on the right includes: Bladder, patient, PTV_4500, PTV_5000, Rectum, Z1-Bridge, Z10-Couch Support, Z2a-Bridge, Z2b-Bridge, Z3-Bridge, Z4-Couch Support, Z5-Hard plate, Z6-Couch Support, Z7-Couch Support, and Z8-Mattress. The Structures table at the bottom lists the following data:

Name	Color	Visible	Volume (cm ³)	Type	Force ED	Fill ED	Relative ED	Show 2D Outlines	2D Transparency	3D Transparency
Bladder	Red	☑	416.647	Internal	☐	☐		☑		☑
patient	Blue	☑	15087.384	External	☐	☐		☑		☑
PTV_4500	Magenta	☑	1144.74	Internal	☐	☐		☑		☑
PTV_5000	Green	☑	72.747	Internal	☐	☐		☑		☑
Rectum	Orange	☑	29.288	Internal	☑	☑		☑		☑
Z1-bridge	Blue	☑	230.411	Couch	☑	☐	1.674	☑		☑
Z10-Couch Support	Blue	☑	272.421	Couch	☑	☐	0.141	☑		☑
Z2a-Bridge	Blue	☑	1501.902	Couch	☑	☐	1.674	☑		☑
Z2b-Bridge	Blue	☑	1119.439	Couch	☑	☐	0.069	☑		☑

The interface also shows a 'Motion Monitoring' panel on the right with buttons for 'Start Motion Monitoring' and 'Stop', and a 'Finish & Continue' button. The bottom status bar displays: Slice Mode S (x): -0.51; T (y): -2.65; C (z): -0.66; cm 82; Isocenter Shift (cm): X: -0.15; Y: -3.14; Z: -0.98; Adaptation Method: Optimize Weights and Shapes.

Unity Applications Software

Clinical Administration: Applying Activity Templates during the Online Workflow

The screenshot displays the Monaco software interface for a treatment session. The top-left panel shows patient information: "Rectum, One", MR#: rectumumb1, 12/12/1970, 47 yr, and 154.1* Rectum, NOS Affirmed. The main workspace is divided into three viewports: a top-left axial view, a top-right sagittal view, and a bottom-left axial view. A red circle highlights the rectum in the axial views, and a green outline indicates the monitored structure. The right-hand panel, titled "Image Acquisition", has tabs for "Pre-Treatment", "Motion Monitoring", and "Post-Treatment". Under "Motion Monitoring", the "Preset" is set to "bFFE Cor Sag Tra RealTime" and the "Monitored Structure" is "Rectum". Buttons for "Start Motion Monitoring" and "Stop" are visible. The bottom panel shows the "Fusion" workflow with a "Reference Plan" of "CT1: 01Rectum" and "Session Images" of "MR&T: No Descrip". A red circle highlights the fusion area in the bottom-left viewport. The bottom-right panel shows a "Structure Visibili" list with various structures and their visibility status.

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The image displays several overlapping screenshots of the MOSAIQ software interface, illustrating the workflow for patient registration, assessment, diagnosis, and radiation prescription.

Patient Registration: Shows fields for Name (Last: ONE, First: CIRS), Middle, Other, and Salutation.

Initial Assessment: Shows MR Initial Assessment - MR#: 4DCIRS1 ONE, CIRS. Includes fields for Record Date and Time (5/14/2018 10:51:10 AM) and Pacemaker.

Diagnosis: Shows Diagnosis Affirmed: 157.9 Pancreas, NOS. Category: PANCREAS, 157. Includes fields for Wait Code, Initial Diagnosis, Consultation (3/6/2018), Inactive Diagnosis, and Tumor Markers.

Course and Intent: Shows Dx: IIB: *Pancreas, NOS Neoplasm, benign. Includes fields for Course (1), Intent (1_Curative), Care Plan (Rad Onc Complex (default)), and Regimen.

Radiation Prescription: Shows Rx Site: Pancreas, Status: Pending, Technique: MR/RT, Modality: Xrays. Includes a table for Dose Specification:

Rx Dose	Fractional Dose	Number of Fractions	Fractionation Pattern	Status
7,400 cGy	200 cGy	37	Daily	

Additional fields include Comment, Close Display Folder, and Care Plan Will Be Changed.

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Patient Care - Plan of Care

ONE, CIRS
MR#: 4DCIRS1
2/2/1992
26 yr

Intent and Diagnosis: Cure - 157.9* Pancreas, NOS Neoplasm, benign, T1, N1, M0, IIB
• Pancreas or hepatopancreatic ampulla

New Intent

Plan of Care

New Activity

Abdomen-HFS-Motion[btFFE Cor Sag Tra RealTime]

Abdomen-HFS-Pre[T2 3D Tra]

Activity Association to Prescription

Activity	Pancreas 7,400.00 cGy / 37# 200.00 cGy Other
Abdomen-HFS-Pre[T2 3D Tra]	<input checked="" type="checkbox"/>


Image Acquisition

Modality: MRI Exam Card: MRL Abdomen HF Acquisition: T2 3D Tra Purpose: Plan Adaptation

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Electron Density Info: Reviewed
Not Reviewed
Reviewed

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MR Daily Assessment - MR#: 4DCIRS1 ONE, CIRS

Info MAR Summary Note Status Print Save Cancel

Record Date and Time: 5/14/2018 12:40:58 PM

Has anything changed?

Pacemaker:

Implants, wires or other:

Metal particles in your eyes:

Shrapnel or bullets in body:

Breast expander w/ magnet:

Hydrocephalus/insulin/med pump:

New Surgery since last visit:

Tattoos or skin patches:


Colored contact lenses:

Pregnant:

Unity Applications Software

MOSAIQ: Treatment Session Manager (Setup and Image)

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




Unity Applications Software

Treatment Session Manager (Motion Monitoring)

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N/A Fields:
Treated Fields:
Image Only:

Unity Applications Software Review Question

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MR Initial Assessment登録

Adapt to Shape or Adapt to Position workflowを実行

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


Unity Applications Software Summary

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1. Clinical Administration
2. Plan of Care
3. MOSAIQ
4. Offline Monaco and Online Monaco
5. Treatment Session Manager(TSM)

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