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
The screenshot shows the AQUA dashboard interface. At the top, there are navigation tabs: Physical, Tests, Schedule, Equipment, System Data, Compliance, and Help. The main header displays the AQUA logo and the institution name. Below the header, there is a 'Dashboard' section with a 'Machine Group' dropdown menu set to 'All Inacs' and an 'Enable Refresh' checkbox. The main content is a table with the following columns: Machine, Machine Status, Mandatory, Compliance, and Parts.

Machine	Machine Status	Mandatory	Compliance	Parts
SA04	Working - Clinical	Yellow	Yellow	Green
EA05	Working - Clinical	White	White	Green
NA00	Working - Clinical	White	White	Green
NA11	Working - Clinical	White	White	Green
WW15	Working - Clinical	White	White	Green

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### Tests

Test List Test Groups File Upload Test Merge

Test List > Edit Test > Edit Test Instance

ID 10472 Name Output Flatness Symmetry - Test VersaHD (6 MV) Equipment/Machine Test VersaHD Complete false

+ Procedure

+ Specification

SAVE

DELETE

CANCEL

### Parameters

Label	Name	Type	Value
Conversion Factor	conversionFactor	test	
Delivery Energy	energy	equipment	6 MV
Flatness Fail Tolerance	flatnessFailTolerance	test	
Flatness Warn Tolerance	flatnessWarnTolerance	test	
Output Baseline	baselineValue	test	
Output Fail Tolerance	outputFailTolerance	test	3.0 %
Output Warn Tolerance	outputWarnTolerance	test	
Radial Flatness Baseline	radialFlatnessBaseline	test	
Radial Symmetry Baseline	radialSymmetryBaseline	test	

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### Tests

- Test List
- Test Groups
- File Upload
- Test Merge

#### Test List > Edit Test

ID	2357	Name	Output Flatness Symmetry	Version	9	Provenance	Demo System - Do not use in a clinical environment	<input type="checkbox"/> Draft
+ Description								
+ Script								
+ Procedure								
+ Specification								
+ Test Images								
Parent Tests	-							

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Test Name	Period	Performance	Compliance Status
Audiovisual Monitors Intercom	Daily	—	DUE
Audiovisual Monitors Video Monitor	Daily	—	DUE
Beam On Indicator Control Warning Lights	Daily	—	DUE
Beam On Indicator Room Warning Lights	Daily	—	DUE
EPID and XVI Touchguards	Daily	—	DUE
Radiation Monitor	Daily	—	DUE
Treatment Door Closing Safety	Daily	—	DUE
Treatment Room Door Interlock	Daily	—	DUE
Accessory Mounting Latches	Monthly	—	DUE
In Room Respiratory Monitoring System	Monthly	—	DUE

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