

Clinic: _____

Date: _____

GAMMEX RMI 467
S/N: 802428-1240

Scanner ID: _____

撮影条件:

ROD#	CT値	rED	Label
		0.01	Air
2		0.258	LN-300 Lung
8		0.424	LN-450 Lung
11		0.926	AP6 Adipose
16		0.954	BR-12 Breast
1		0.988	CT Solid Water
3		0.988	CT Solid Water
5		0.988	CT Solid Water
13		0.988	CT Solid Water
15		1.000	Water Insert (容器)
10		1.049	BRN-SR2 Brain
6		1.065	LV1 Liver
4		1.098	IB Inner Bone
7		1.106	B200 Bone Mineral
14		1.280	CB2-30% CaCO3a
12		1.470	CB2-50% CaCO3a
9		1.697	SB3 Cortical Bone

並び変えるときは自分で番号を記入のこと

