

Clinic: _____

Date: _____

GAMMEX RMI 467
S/N: 802428-1124

Scanner ID: _____

撮影条件:

ROD#	CT値	rED	Label
		0.01	Air
2		0.267	LN-300 Lung
8		0.416	LN-450 Lung
11		0.930	AP6 Adipose
16		0.957	BR-12 Breast
1		0.988	CT Solid Water
3		0.988	CT Solid Water
5		0.988	CT Solid Water
13		0.988	CT Solid Water
15		1.000	Water Insert (容器)
10		1.047	BRN-SR2 Brain
6		1.072	LV1 Liver
4		1.097	IB Inner Bone
7		1.105	B200 Bone Mineral
14		1.273	CB2-30% CaCO3a
12		1.466	CB2-50% CaCO3a
9		1.695	SB3 Cortical Bone

並び変えるときは自分で番号を記入のこと

